



NeoLumo[®]
C O L L E C T I O N

ObeyYourBody[®]
NeoLumo[®]
C O L L E C T I O N



neolumo.com



NeoLumo
COLLECTION

Contents

Everything you need to know about the miraculous world of NeoLumo®

Welcome to the Amazing world of NeoLumo®	1-2
The Innovative NeoLumo® Collection	
Our LED Light Technology	3-4
Our Other Technologies	5-6
Skin Cells Regeneration Changes Over Time	7-8
Our Products	9-10
NeoLumo® Integrated Skincare Systems	
Clarity Forehead Treatment	11-12
Renewal Cheek & Around Lip Treatment	13-14
Balance Neck Treatment	15-16
Shaping Eye Treatment	17-18
How to use Neolumo® Integrated Systems according to skin's needs	19-20
NeoLumo® Skincare Electronic Devices.....	21-22
Facial Multi-Therapy Device	23-26
Facial Infrared Renewal Device	27-30
Neck Lift Premium Device	31-32
Neck LED Spectrum Device	33-34
Caviar Royal Mask	35-36
Recommended by Medical Professionals / Limited Lifetime Warranty.....	37-38
Case Studies Results Say It All	39-40

NeoLumo®
C O L L E C T I O N



Welcome to the Amazing world of NeoLumo®

ObeyYourBody® is a globally acclaimed health and beauty cosmeceuticals brand with over 15 years of successful international track record. Our creams and serums are formulated with the best that nature and science have to offer.

In the last few years a new and very promising field in skincare has been developed, indirectly originated by no other than NASA: when astronauts go out into space, their cells stop multiplying because there is no gravity, sunlight, or oxygen - so if they got a cut or a bruise, it wouldn't heal. They needed something that would help their cell regeneration, so if anything happened, their body could heal. Thus, LED light therapy came into being. This breakthrough technology actually was originally initiated by NASA for plant growth experiments on shuttle missions.

So, LED was created for astronauts and then the health industry got a hold of it and started to develop electronic devices that use LED lights in different wavelengths and carry profound benefits to the human skin.

In some cases, LED lights have been combined with other beneficial technologies such as RF (Radio Frequency) and EMS (Electronic Muscle Stimulation) and more.

These advanced devices are highly effective in combating aging skin, wrinkles, blemishes, acne, pigmentation and other skin concerns. Another benefit is quicker and deeper absorbance of creams and serums applied to the skin. The NeoLumo® devices are safe to use and have certificates depending on the specific device: FDA, CE, RoHS, ISO 9001.

ObeyYourBody® is on the forefront of this skincare revolution, developing and bringing to you, the consumer, the absolute best skincare devices on the market. And we do not stop there. Due to our heritage, we are situated in a unique position to combine the amazing skincare benefits of these electronic devices with specially formulated creams and serums to offer perfect synergy of both and allow for beautiful, young looking and refreshed skin.

In this book we've shared more details about this new skincare revolution: technology, synergies, results and more interesting facts. We promise to continue to be your ultimate skincare brand, relentlessly bringing to you the most effective skincare solutions.

Join our countless satisfied customers worldwide and enjoy your new skincare regime!

The Innovative NeoLumo® C O L L E C T I O N

ObeyYourBody®



Uniquely targets specific facial zones

Light Therapy | Active Cream | Vibration Osmosis

Renewal Cheek & Around Lip Treatment

Strengthens & Smooths
GREEN Triple-Action
Unique Anti-Pigmentation System



Clarity Forehead Treatment

Boosts Collagen
RED Triple-Action
Unique Anti-Aging System



Balance Neck Treatment

Firms & Tightens | BLUE Triple-Action
Unique Cleansing & Antibacterial System




Shaping Eye Treatment

Fades Fine Lines & Dark Circles
Innovative Technology
Unique Dual-Action Eye Contour Treatment



Facial Multi-Therapy Device

5-in-1 active anti-aging system



- Stimulates collagen
- Smooths lines & wrinkles
- Treats blemishes & UV damage

Neck LED Spectrum Device

Triple lift & revitalizing care



- Stimulates collagen
- Smooths & balances
- Firms up sagging skin

Neck Lift Premium Device

Triple lift & revitalizing care



- Lifts & tightens skin
- Reduces lines & wrinkles
- Promotes nutrient absorption

Facial Infrared Renewal Device


Anti-aging cell rejuvenation



- Boosts natural healing
- Stimulates cell renewal
- Reduces pore size

Skin Cells Regeneration Changes Over Time

COLLAGEN IN HUMAN SKIN




The human skin is the largest body organ. The top part and the thinnest is the **Epidermis**. The upper part of the Epidermis consists of dead cells that are consistently removed and replaced by new cells. This process slows down with age and environmental radiation and prolusion. The Epidermis is responsible for our skin tone due to the concentration of Melanocytes. The Epidermis has also an important part in protection against microorganisms.


Underneath the epidermis lies the **Dermis**. The dermis contains blood vessels that nourish the skin with oxygen and nutrients. The blood vessels also allow immune system cells to aid the skin to fight infections. These vessels also help carry away waste products. With aging, these capacities are reduced, causing skin to look more dull and prone to blemishes. The dermis includes also various glands, hair follicles and Collagen, a protein that is responsible for giving skin strength and elasticity. With aging, collagen diminishes and that is the main cause for aged skin, wrinkles, fine lines and sagging.

Beneath the dermis is the deepest layer of our skin, **Hypodermis**. It contains many collagen cells as well as fat. With aging, the volume of facial fat decreases and there is less supportive tissue to support the normal turgor and elasticity of the skin. The facial skin begins to droop and sag resulting in a look that can be interpreted as appearing old and tired. The bones and muscles of the face also lose volume.

Before NeoLumo®



After NeoLumo®



NeoLumo® MAKES ALL THE DIFFERENCE!

Our LED Therapy and other beneficial technologies deeply penetrate the skin: allowing for the generation of collagen and new cells. Wrinkles and fine lines diminish, skin tone dramatically improves, acne conditions greatly improve, and the skin becomes much more elastic, smooth and healthy.

Some of the devices' functions also enhance absorption of our specially formulated creams and serums into the skin. Both working synergistically to provide amazing skincare results!

NeOLumo[®]
C O L L E C T I O N

Our Products

Integrated Skincare Systems

Uniquely targets specific facial zones
Light Therapy | Active Cream | Vibration Osmosis



Unique Triple-Action systems pairing highly-beneficial photon light wave and vibration technology with synergistically matching creams.

Each system is specially formulated to treat a certain common skin concern: wrinkles, hyperpigmentation, sagginess & blemishes, eye area fine lines & wrinkles.

Each system is fully suitable to be used on all areas of face and neck, however each system has been designated to a specific facial zone.

Each system comes in a fully integrated single compact unit!

Clarity Forehead Treatment

Boosts Collagen | **RED Triple-Action**
A Unique Anti-Aging System



A highly effective Anti-Aging treatment all in one amazing product!

To treat wrinkles and fine line on all areas of face and neck. Especially recommended for the forehead's deep wrinkles.


Synergetic treatment of red LED, vibration and active ingredients that visibly reduce lines and wrinkles.

Stimulates collagen and elastin production, increases cells' energy levels, hydrates, and repairs sun damage.

The red LED works in synergy with the active ingredients in the cream, while the vibration function helps the skin to absorb the active ingredients, altogether enhancing collagen production, accelerating cell renewal and creating an anti-aging effect that reduces blemishes and unifies the skin tone for firm and young look.


RED LED LIGHT

- Anti-aging - helps to reduce wrinkles.
- Stimulates collagen and elastin production, reducing fine lines and wrinkles.
- Increases blood circulation, helping absorption of cream's active nutrients into the skin's layers.
- Increases production of ATP (energy engine of cells).
- Antioxidant, helps repair sun damage.
- Helps lighten and diminish scars and even overall skin tone.




VIBRATION MASSAGE

- Vibrates 7000 times per minute.
- Boosts absorption of the cream's active ingredients.
- Accelerates blood circulation, helps to produce collagen and elastin.
- Reduces lines and wrinkles.
- Aids in lymphatic drainage.
- Cleansing and exfoliation effect, ousting dirt and debris.



ACTIVE INGREDIENTS OF THE CREAM

- Syn-Coll and Matrixyl® 3000 peptides encourage collagen formation.
- Moisturizing Marula Oil and Hyaluronic Acid support an optimal hydration balance.
- Vitamin C promotes a brighter complexion.
- Enriched with antioxidant Vitamin E to combat harmful free radicals.



"I love this Neolumo product. I started to see results from the very first time I used it, and it just gets better...This is magical!"

Suzanne E.

Renewal Cheek & Around Lip Treatment

Strengthens & Smooths | **GREEN Triple-Action**
A Unique Anti-Pigmentation System



A highly effective anti-pigmentation treatment all in one amazing product!

To treat hyperpigmentation on all areas of face and neck.
Especially recommended for cheeks and around lips areas.


Synergetic treatment of green LED, vibration and active ingredients that balance melanin in the skin.

Fades pigmentation marks, balances and blends contours to a natural and uniform skin tone.

The green LED works in synergy with the active ingredients in the cream, while the vibration function helps the skin to absorb the active ingredients, and altogether they balance and catalyze cell renewal and metabolism and reduce melanocyte activity, for smooth, young and even tone skin.


GREEN LED LIGHT

- Balances, calms and renews.
- Increases cellular renewal and energy.
- Calms damaged skin.
- Inhibits the production of excess melanin, prevents hyper-pigmentation - for even skin tone.
- Stimulates absorption of cream's active ingredients into the skin layers.



VIBRATION MASSAGE

- Vibrates 7000 times per minute.
- Boosts absorption of the cream's active ingredients.
- Accelerates blood circulation, helps to produce collagen and elastin.
- Reduces lines and wrinkles.
- Aids in lymphatic drainage.
- Cleansing and exfoliation effect, ousting dirt and debris.



ACTIVE INGREDIENTS OF THE CREAM

- Matrixyl® Morphomics and SWT-7™ smooth the appearance of smile lines.
- Moisturizing Marula Oil and Hyaluronic Acid support an optimal hydration balance.
- Essential ATP restores energy to the skin, for an overall anti-aging effect.
- Enriched with antioxidant Vitamin E and skin-smoothing Vitamin B5.



"My marionette and "barcode" lines decreased and the area around my lips and cheeks looks so much better. I like now to look at myself in the mirror"
Betty W.

Balance Neck Treatment

Firms & Tightens | BLUE Triple-Action
A Unique Cleansing & Antibacterial System



A highly effective Cleansing & Antibacterial treatment - all in one amazing product!

To treat Sagginess and blemishes on all areas of face and neck.
Especially recommended for the neck.

Synergetic treatment of blue LED, vibration and active ingredients that disinfect the skin.

Cleanses and destroys the bacteria present in all layers of the skin, including acne causing bacteria.

Destroys black and white spots and prevents their reappearance on the skin.

The blue LED works in synergy with the active ingredients in the cream, while the vibration function helps the skin to absorb the active ingredients, altogether encouraging cell growth along with antibacterial activity to treat blemishes, for a smooth and supple skin with an amazing pushup effect.

BLUE LED LIGHT




- Soothes, calms and tightens loosened skin.
- Reduces excess oil production.
- Tightens loosened skin.
- kills acne-causing bacteria beneath the skin, to prevent and treat acne.
- Decreases rashes and scars.

VIBRATION MASSAGE

- Vibrates 7000 times per minute.
- Boosts absorption of the cream's active ingredients.
- Accelerates blood circulation, helps to produce collagen and elastin.
- Reduces lines and wrinkles.
- Aids in lymphatic drainage.
- Cleansing and exfoliation effect, ousting dirt and debris.

ACTIVE INGREDIENTS OF THE CREAM

- Rich in Zinc, Copper & Magnesium and Septonic M4 that reconstruct the skin barrier.
- Moisturizing Marula oil and hyaluronic acid support an optimal hydration balance.
- Innovative Voluform™ offers a remarkable skin firming and plumping effect.
- Enriched with Vitamin A and antioxidant Vitamin E to enhance skin radiance.



"Ever since I started to use the NeoLumo product, I enjoy smooth and clear skin clear of blemishes and breakouts."

Dianne K.

Shaping Eye Treatment

Fades Fine Lines & Dark Circles | Innovative Dual-Action

A Unique Complete Eye Contour System



**A highly effective eye area shaping treatment:
All in one amazing product!**

To treat the Eye Area: Crow-feet and under-eye bags and dark circles.

A special Dual-Action treatment system designed for the delicate eye contour and under eye area.

The integrated vibration function gently soothes fine lines, reduces wrinkles, and increases the skin's absorption of the cream's active ingredients.


The cream's active ingredients, work in synergy with the vibration function for improved absorption, strengthen the fragile microvascular infrastructure, raise the upper eyelid, and reduce "crow's feet", for a youthful looking smooth skin.

VIBRATION MASSAGE

- Vibrates 7000 times per minute.
- Boosts absorption of the cream's active ingredients.
- Accelerates blood circulation, helps to produce collagen and elastin.
- Reduces lines and wrinkles.
- Aids in lymphatic drainage.
- Cleansing and exfoliation effect, ousting dirt and debris.

ACTIVE INGREDIENTS OF THE CREAM


- Beautifye™ helps to lift the upper eyelid, reduce crow's feet, and diminish puffiness. Also reduces pigments responsible for uneven coloration and dark circles.
- The Syn-Coll and Matrixyl® 3000 peptides encourage collagen formation.
- Moisturizing Marula Oil and Hyaluronic Acid support an optimal hydration balance.



Beautifye™ is an award-winning ingredient manufactured by the French company Sederma that directly targets the skin of the upper eye, supporting healthier skin that's less likely to sag or show wrinkles.

105 female volunteers, mean age 51 (40 to 60 years old), with a tired look characterised by wrinkles, puffiness and dark circles, applied twice daily for 28 days a cream containing 3% BEAUTIFEYE™ in the area around the eyes. Consumer test conducted by Eurofins.

Within 7 days only, BEAUTIFEYE™ influences the way volunteers perceive themselves. The aspect of fatigue of the eye area is significantly reduced.



"I'm so happy with this Neolumo product! The fine lines and dark circles are almost gone now. I'll keep using this product and highly recommend it to everyone concerned about their eye contour"

Greta M.

How to use NeoLumo® Integrated systems according to skin's needs

Anti-Aging



Anti-Pigmentation + Anti-Aging



Anti-Acne + Anti-Aging



All-Round Anti-Aging + Anti-Pigmentation + Anti-Acne



+



+



+



+




NeoLumo[®]

COLLECTION

Electronic Skincare Devices

Facial Multi-Therapy Device

5-in-1 active anti-aging system



Neck LED Spectrum Device

Triple lift & revitalizing care



Neck Lift Premium Device

Triple lift & revitalizing care



Facial Infrared Renewal Device


Anti-aging cell rejuvenation



Renew your youth with the gentle smoothing qualities of natural light waves coupled with additional proven effective technologies. This breakthrough concept from ObeyYourBody twins the eternal wonders of nature with state-of-the-art skincare science. Originally developed by NASA scientists, this anti-aging tech innovation is now being applied to treat common signs of aging and skin disorders. NeoLumo presents a genuine, comprehensive skin transformation, for a firmer, smoother complexion with a dramatically enhanced texture and a clearer, more youthful appearance.



Facial Multi-Therapy Device 5-in-1 active anti-aging system



This state-of-the-art device leverages 5 essential skincare concepts to present the ultimate in effective, non-invasive skin rejuvenation:

Electrical Muscle Stimulation (EMS) to stimulate collagen and reduce wrinkles.

Mesoporation to deliver active ingredients to every layer of the skin.

Electroporation to transport hyaluronic acid and collagen even deeper into the skin.

Radio Frequency (RF) to boost collagen production.

5 LED lights therapy to treat wrinkles, blemishes, pigmentation and sun damage.

Together, these breakthrough technologies gently yet thoroughly cultivate a fresher, more radiant, and beautifully clarified complexion with a smooth, unified tone that shines with renewed youth.



KEY BENEFITS:

- Stimulates collagen
- Smooths lines & wrinkles
- Imparts vital nourishment
- Treats blemishes & UV damage
- Visible anti-aging effect



"I purchased a Neolumo multi therapy device and was quite skeptical. Look at me now, what a surprise! My skin looks smooth, clear, even, tight, and seems much younger.

Thanks Neolumo!"

Vanessa J.

5 highly beneficial technologies in one device

1. LED Lights

Originally developed by NASA for astronaut's cell regeneration and growing plants in space, and then used by the health industry with a wide range of amazing benefits.




Red LED light therapy is used to help regenerate skin cells.

Blue LED is used for Acne and blemishes.


Green LED is used for pigmentation control.

Yellow and **Pink** LEDs are used for cleansing and purifying the skin.




3. Mesoporation

Mesoporation comes from the combination of mesotherapy and electroporation. Mesotherapy is a system of transdermal skin delivery in which micro-needles puncture the skin to allow the delivery of active ingredients to deeper layers of the skin. Electroporation is a non-invasive electrical treatment which also carries active ingredients into the dermal layers. Unlike mesotherapy, mesoporation does not require needles.




4. Electroporation

Electroporation provides a means to transport collagen and hyaluronic acid deeper into the skin. Electroporation provides the means to achieve better results with your skin care products without the need for needles.




2. Ems (electrical muscle stimulation)

Ems is used to activate the skin cells, stimulating rejuvenation of collagen, to improve skin elasticity, reduce wrinkles and leave skin looking smooth and healthy.



5. Radio Frequency (RF)

Radio Frequency energy uses the tissue's natural resistance within the various layers of the skin to transform the RF energy applied to the skin into thermal energy. RF energies can be used for all skin types and in different depths of penetration, allowing for increased collagen contraction and production of new collagen.



Facial Infrared Renewal Device

Anti-aging cell rejuvenation



Red Infrared LED light therapy represents a genuine skincare revolution, transforming daily facial regimens and far surpassing expectations.

This unique device leverages the true power of infrared LEDs to deeply penetrate the skin, converting light energy into cellular energy and rapidly stimulating natural healing processes. As new skin cells form, wrinkles and fine lines are gently smoothed away, and unsightly scars and pigmentation spots are replaced by a clear, unified complexion.

Without any need for painful, invasive procedures, the face is effortlessly beautified, presenting a glowing, youthful look that is significantly firmer, smoother, and radiantly flawless.



"I was concerned about my wrinkles, pigmentation and age spots. Only a few weeks using the NeoLumo Infrared device, all of these got so much better. I'm delighted and very happy with the results and will continue to use it."


Judy A.

KEY BENEFITS:

- Boosts natural healing
- Increases collagen production
- Stimulates cell renewal
- Reduces pore size
- Beautifies the appearance
- Removes wrinkles and fine lines
- Visible anti-aging effect


INTERNAL SKIN TREATMENT

The infrared light penetrates the skin and helps cells regenerate or repair themselves. The surface might feel warm.




SUPPORTS SKIN REJUVENATION

Quickly stimulates natural healing of the skin.




REBUILDS SKIN CELLS

New skin cells will be formed, as wrinkle lines and dark spots disappear.



The device heats up to 42° Celsius, to increase the blood flow to the treated area, which helps the LED lights to penetrate even deeper and quicker into your skin.




Originally developed by NASA for astronaut's cell regeneration and growing plants in space, and then used by the health industry with a wide range of amazing benefits.




The science of LED lights


different depths of penetration



Wavelength (nm)	Color Range	Penetration (mm)
150-380	UV	<0.1
390-470	Violet to Deep Blue	~0.3
475-545	Blue-Green	~0.3-0.5
545-600	Yellow-Orange	~0.5-1.0
600-650	Red	~1.0-2.0
650-950	Deep Red-Near Infrared	2-3
950-1200	Infrared	1




SKIN BEFORE



PERFECTIO AT WORK

Facial Infrared Renewal Device – at Work
Deep skin penetration treatment




SKIN AFTER

- Strong light exposure
- Gentle and painless
- Smooth operation and application

Neck Lift Premium Device

Triple lift & revitalizing care



Anti-Aging - Collagen Boosting - Lead-in Nutrition Waterproof - Automatic Operation

Featuring an advanced aesthetic design and multiple innovative technologies, this anti-aging neck therapy device presents the ideal solution for common issues like wrinkles, a double chin, and sagging skin. The impressive micro-arc massage head is designed to fit perfectly and deliver superior care in three special modes. The **Normal** mode vibrates and massages to enhance circulation and skin vitality. The **Cool** mode uses blue light to lift and tighten the skin while relieving muscle fatigue and dark circles. Finally, the **Hot** mode draws on the powers of red light to reduce wrinkles, boost radiance, and promote nutrient absorption for a visibly revitalized complexion. Truly the optimal complete solution for renewed youthful firmness and a flawlessly sculpted look.





"The Neolumo neck lift device definitely kept its promise. My neck looks so much better now, lifted, firm, clear and much younger."
Bella Z.

Red LED light therapy is used to help regenerate skin cells.
Blue LED is used for acne and blemishes.



KEY BENEFITS:

- Lifts & tightens skin
- Reduces lines & wrinkles
- Enhances youthful radiance
- Promotes nutrients absorption



Neck LED Spectrum Device Dynamic collagen boosting



Ergonomically designed to pamper the neck and shoulder areas, this outstanding therapy device features three advanced light massage modes to imbue every part of the skin with rejuvenating energy.

Blue LED fights blemishes and acne. **Green** LED light smooths and balances the skin, reducing redness while enhancing cell metabolism. **Red** LED light boosts cell vitality and circulation, minimizing the look of wrinkles, scars, and pigmentation spots.

The device is also ideal for an acupuncture-inspired neck massage, helping to relieve stress and muscle tension while simultaneously promoting a youthful glow. It's a full spectrum of revitalizing care, ensuring you look – and feel – rested, refreshed, and radiantly beautiful.




"My NeoLumo neck LED device is all I need for my neck. I love it!"
Julie R.

Red LED light therapy is used to help regenerate skin cells.
Blue LED is used for acne and blemishes.
Green LED is used for pigmentation control.



KEY BENEFITS:

- Stimulates collagen
- Smooths & balances
- Firms up sagging skin
- Promotes nutrient absorption
- Visible anti-aging effect



Caviar Royal Mask
Age-Reverse Collagen Treatment



Anti-Aging - Collagen Boosting - Natural UV Protection
World-Class Ingredients - Perfect Synergy with
NeoLumo Devices

Give your skin the optimal age-reverse treatment with our Caviar Royal Mask. This unique mask contains two of the most powerful age-reverse agents: **caviar and collagen**, complemented with rare minerals and herbal extracts. The Caviar Royal Mask will make your skin firmer, smoother, and rejuvenated with luminous glow, leaving your face looking instantly younger and radiant.

Caviar Extract benefits to the skin

When you think of Caviar, the first connotation is expensive delicacy of the rich and famous. But caviar actually has highly beneficial benefits for the skin. It is a powerful anti-aging agent. The caviar in our mask has strong antioxidant properties which help prevent collagen and elastin breakdown, while eliminating wrinkles and fine lines.

Caviar is also a strong moisturizer, containing fatty acids Omega 3 and Omega 6, which have an important role in maintaining the integrity and barrier function of the skin. It is rich in vitamins, proteins, nucleic acids and phospholipids, to make the skin elastic, firm, with velvety texture.


It contains potassium, sodium and calcium that enhance the skin's luminosity and feed the cells. Caviar extract is ideal to repair skin damage as well as stimulate the metabolism of the skin cells.



Collagen benefits to the skin

Collagen is the most abundant protein in the skin and makes up 75-80% of your skin. Collagen gives the skin its elasticity and strength, contributing to a youthful appearance.

That's why it is so important. However, with aging, your body's collagen levels go down. The Collagen in our mask helps compensate for this collagen loss, working on the surface of the skin while also stimulating collagen production in its deeper layers. It enhances hydration by reducing water loss, and also adds a barrier against environmental damages. It helps remove scars and wrinkles, while helping your skin cells renewing and repairing, making your face smooth, glowing, firm, and younger looking.



USE IN PERFECT SYNERGY WITH THE NEOLUMO [®] LED LIGHTS TO ENHANCE EFFECTIVENESS		
LED Light	LED Main Benefit	Mask Synergetic Ingredients
Red	Anti-Aging	Collagen, Caviar Extract, Green Tea Extract
Green	Anti-Pigmentation	Licorice Extract, Baikal Skullcap
Blue	Calming and Healing	Rosemary Extract, German Chamomile, Oat Extract
Yellow + Pink	Cleansing and	Sweet Almond Oil, Japanese Knotweed, Kaolin

Recommended by Medical Professionals

Dr. Alex Levin

"Research confirms that LED lights and other technologies such as RF and EMS are highly effective in combating the signs of aging skin. They increase collagen formation, reduce wrinkles and repair damaged skin. They also treat pigmented skin conditions and treat pores for a more balanced skin tone. Another benefit is the betterment of acne prone skin, acting against blemishes and brightening the skin.

Using these devices in conjunction with the specially formulated creams and masks further improves the visible results for a healthier, smoother and firmer complexion."



NeoLumo®

C O L L E C T I O N

NeoLumo's high quality products meet strict international standards and regulations. This is why we can **provide you with a Limited Lifetime Warranty for our wide range of devices.**



ObeyYourBody®

NeoLumo
C O L L E C T I O N

Limited Lifetime Warranty

We offer limited lifetime warranty for your NeoLumo® Device. If your product has manufacturing or workmanship defects, ObeyYourBody® will repair or replace it. Warranty does not cover normal wear and tear of the device. Warranty does not cover misuse or mistreatment of the device. You must register the device in 7 business days after purchase date at www.obeyyourbody.com in order for the warranty to be valid.



Case Studies Of NeoLumo[®] Users Say It All

96%

agree their wrinkles and fine lines reduced/diminished.

100%

agree their skin looks fresher and younger.

98%

agree their skin looks visibly firmer and smoother.

97%

agree their complexion is more clear and even.

93%

agree that their acne condition significantly improved.

95%

agree that their skin pigmentation condition has been greatly improved.

100%

agree that NeoLumo is the most effective skincare solution they ever used.

100%

agree they will recommend NeoLumo to others.